


How Does The PractiVet Infusion Pump Compare to Other Leading Veterinary Pumps?



Vet Pro 2000



Heska



PractiVet

Price	~\$1100	~\$1200	\$895
Warranty	12 Months	12 Months	24 Months*
Loaner Pump Program	No	No	Yes*
Dimensions	6.75x3.5x3.5 in	5.4x5.4x2.8 in	5.7x4.7x3.9 in
Weight	1 lb	2.6 lbs	2.2 lbs
Battery	1800 mAh Li-Ion	Nickel-Cadmium 0.9 Ah	1900 mAh Li- Polymer
Display	LCD	LCD	Color LCD
Programmed Admin Sets	1	1	10
Infusion History	Last 24 Hours	None	1300 Records
Recalibration	Must Send In Every 1500 Hours Cost ~\$250	Must Send In Every 1500 Hours Cost ~\$250	In Clinic In Less Than 10 Minutes
Double CPU	No	No	Yes
Ability To Upload Data	No	No	Yes
Wifi Connectivity	No	No	Yes

*When using PractiVet fluid therapy devices, you gain an additional 12 month warranty and access to our loaner pump program