


FASTFLASH

THE NEW RANGE OF IV CATHETERS WITH EARLY FLASH BACK

- It provides an early sign of vein reach thanks to a notch on the distal part
- Immediately gives a feed back about blood flow
- Flash back confirmed by the chamber



① **FIRST FLASH**

NOTCHED NEEDLE ALLOWS EASY, IMMEDIATE FLASHBACK

③ **SECOND FLASH**

FLASHBACK CHAMBER

②

THE TRANSPARENT VENOUS CATHETER ALLOWS AN IMMEDIATE CLEAR VIEW



"Partnering with you to bring New Technology, Quality, Reliability and Service Above and Beyond"