




### Standard IV Catheters

**Features:**

- Very economical IV catheter
- Triple bevel needle for smooth insertion
  - Bevels form a blade which makes a small, surgical incision
  - Reduces force needed to penetrate the skin
  - Also reduces pain and tissue trauma
- Catheter is pre-lubricated to more easily slide into the vein
- Designed for extended intravenous use
- Clear flashback chamber
- Luer Lock cap for the closing of the proximal end of the catheter after removal of the needle



SKU#	DESCRIPTION	PACK
059162	Standard IV Catheter, Orange, 14 G x 2 In	EACH
059171	Standard IV Catheter, Green, 18 G x 1.25 In	EACH
059172	Standard IV Catheter, Green, 18 G x 2 In	EACH
059173	Standard IV Catheter, Pink, 20 G x 1 In	EACH
059174	Standard IV Catheter, Pink, 20 G x 1.25 In	EACH
059175	Standard IV Catheter, Blue, 22 G x 1 In	EACH
059176	Standard IV Catheter, Yellow, 24 G x 0.75 In	EACH



### Winged Infusion Sets

Winged infusion sets, also called scalp vein or butterfly sets, are designed for IV injection or blood collection.

**Features:**

- Sharp needle for less pain
- Protected luer lock connection
- Economical

SKU#	DESCRIPTION	PACK
059177	Winged Infusion Set, Butterfly Style, Yellow, 19 G x 0.75 In	EACH
059178	Winged Infusion Set, Butterfly Style, Green, 21 G x 0.75 In	EACH
059179	Winged Infusion Set, Butterfly Style, Black, 22 G x 0.75 In	EACH
059180	Winged Infusion Set, Butterfly Style, Blue, 23 G x 0.75 In	EACH
059181	Winged Infusion Set, Butterfly Style, Orange, 25 G x 0.75 In	EACH